

## Need to Replace Your Aging IL943 or Eppendorf ELEX 6361 Flame Photometer?

### Available Services:

- Installation and Training
- IQ/OQ/PQ Services
- Application & Method Development Assistance
- Preventative Maintenance Plans (PMs)
- SOP Writing Assistance
- Let Us Assist You in the Transfer of Your Flame Photometer Applications to the Krüss Flame Photometer Platform



3  
**Li**  
Lithium  
Alkali Metal

11  
**Na**  
Sodium  
Alkali Metal

19  
**K**  
Potassium  
Alkali Metal



The A. Krüss Optronic Flame Photometer is an Ideal Replacement for These Older Flame Photometer Models



### Key Features

- Reliable, Simultaneous Measurement of Na, K, Li and Ca
- Complies with 21 CFR Part 11 and GMP/GLP
- Fast Measurement—Up to 300 measurements per hour
- High Safety Standards Through Intelligent Operational Safety Systems
- Built-in Touchscreen Computer with Easy-to-Use Intuitive Software
- User Administration Feature with Two Authorized Levels
- Four Different Models to Choose From

**Contact Us Today for a Consultation!**