

AGS Scientific Is Now Carrying the Complete Line of A. Krüss Optronic Gemology Products!

Gemology, the science of gemstones and ornamental stones, is a branch of mineralogy. Mineralogists study natural stones as well as artificial gemstones. There are a variety of instruments used by the gemologist to determine the true identification of a stone as well as its authenticity and value.



Key Products

- Spectroscopes and cold light sources to help identify colored stones by their specific absorption spectra
- Refractometers to determine the visual characteristics of ornamental stones and gemstones
- Polariscopes to identify gemstone imitations from real gemstones
- Dichroscopes for the observation of pleochroism in transparent gemstones
- Loupes, tweezers, color filters and UV lamps
- Microscopes for determining the quality of stones as well as for displaying gems and stones to customers
- Diamond Testing instruments

Contact Us Today For More Information!