

## New Camera and Display for Use With A. Krüss Optronic Microscopes and Other Microscope Brands

The MKTV5 Camera makes it easy to connect to a monitor, TV, tablet or smartphone. The standardized USB 2.0 Hi-Speed interface and C-Mount thread makes it easy to connect the MKTV5 Camera to all commercially available microscopes. A monitor is also available that directly mounts on the camera.



### MKTV5 Camera

- 5.0 Megapixel Camera with 1/1.8" CMOS Sensor
- Resolution of 2592 x 1944 Pixels
- USB 2.0 Hi-Speed Interface
- HDMI Interface (HDMI Cable Included)
- Live Image, Image Recording and Video Function
- Storage of Image and Video on SD Card (Included)

### MKTV5 Display

- 11.8" 1080p IPS LC-Display
- Resolution of 1920 x 1080 Pixels
- HDMI Interface for Connection to MKTV5

